

## vita-life® MRT-System

### Contraindications

With regard to contraindications, we distinguish between **absolute contraindications** (in which case the MRT-system may not be used under any circumstances) and **relative contraindications**; if these are present, the MRT-system may only be used under close medical supervision, both for the therapy of the condition itself and for general health care.

#### **ABSOLUTE Contraindications are**

- during pregnancy,
- Children under the age of 2 weeks
- if the person has had electronic device implanted such as pacemaker or insulin pump
- in severe cardiac arrhythmias (from Lown IV b)
- Epilepsie,
- in acute episodes of hyperthyroidism.

#### **RELATIVE contraindications - Use only under medical supervision!**

- Tumoren disease